

# Memo

To: AML Reference Laboratory Clients  
From: Allina Medical Laboratories, Reference Laboratory Services  
Date: October 6, 2008  
Re: Reference Changes, Glomerular Basement Membrane Ab, IgG  
(1016/AGB)

---

Allina Medical Laboratories has received notice from Mayo Medical Laboratories (MML) that **effective November 18, 2008**, the **Glomerular Basement Membrane Ab, IgG (1016/AGB)** will switch from the current ELISA methodology to a Multiplex flow immunoassay method. This conversion will have effects on the reference values for this test

**Reference Values:**

<1.0 U	Negative
≥1.0 U	Positive